

Anti-Sun2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Uh99

Item Number: CM0019108



Introduction

SUN2 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Targets the LINC complex component SUN2, essential for nuclear positioning and migration. Ideal for cytoskeleton and nuclear envelope research.

[Learn More](#)

Parameter	Value
Product Name	SUN2 Rabbit pAb
Remarks/Alias	UNC84B; rab51P; SUN2
Species	Human
GeneID (Human)	25777
GeneID	25777
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 255-572 of human SUN2 (NP_001186508)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UH99
Protein Weight	80kDa
Shipping	Ice bag