


Anti-Suox Rabbit Monoclonal Antibody For Wb And Elisa - P51687

Item Number: CM0003148



Introduction

High-quality rabbit monoclonal antibody against human SUOX (sulfite oxidase). Validated for Western blot and ELISA with cross-reactivity to rat. Ideal for studying sulfur metabolism, mitochondrial function, and hypoxia-related nitrite reduction.

[Learn More](#)

Parameter	Value
Product Name	SUOX Rabbit mAb
Remarks/Alias	SUOX; sulfite oxidase
Species	Human
Gene ID	6821
Immunogen	Synthetic peptide (aa 100-200 of human SUOX, P51687)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P51687
Protein Weight	60kDa
Shipping	Ice bag