


Anti-Supt20H Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Nem7

Item Number: CM0029553



Introduction

Rabbit polyclonal antibody targeting human SUPT20H (SPT20 homolog, p38IP). Validated for Western blot, IHC-P, and ELISA in human, mouse, rat. Applications include autophagy, gastrulation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	SUPT20H Rabbit pAb
Remarks/Alias	C13; FP757; P38IP; SPT20; FAM48A; C13orf19; SUPT20H
Species	Human
GeneID (human)	55578
GeneID	55578
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 580-779 of human SUPT20H (NP_001014308.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NEM7
Protein Weight	86kDa
Shipping	Ice bag