


Anti-Supt5H/Spt5 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - O00267

Item Number: CM0003632



Introduction

Rabbit monoclonal antibody targeting human SUPT5H/SPT5, suitable for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Ideal for transcription elongation and HIV Tat research.

[Learn More](#)

Parameter	Value
Product Name	SUPT5H/SPT5 Rabbit mAb
Remarks/Alias	SPT5; SPT5H; Tat-CT1; SUPT5H/SPT5
Species	Human
Gene ID	6829
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 800-900 of human SUPT5H/SPT5 (O00267).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00267
Protein Weight	121kDa
Shipping	Ice bag