


Anti-Supt6H/Spt6 Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q7Kz85

Item Number: CM0023917



Introduction

High-quality rabbit polyclonal antibody against human SUPT6H (SPT6), validated for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Ideal for transcription elongation and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	SUPT6H Rabbit pAb
Remarks/Alias	SPT6; SPT6H; emb-5; SUPT6H
Species	Human
GeneID (Human)	6830
GeneID	6830
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1500-1726 of human SUPT6H (NP_003161.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7KZ85
Protein Weight	199kDa
Shipping	Ice bag