


## Anti-Supt7L Rabbit Polyclonal Antibody For Wb And Elisa - O94864

Item Number: CM0020277



### Introduction

High-quality rabbit polyclonal antibody against human SUPT7L (STAGA complex 65 subunit gamma), validated for Western blot and ELISA. Reacts with human, mouse, and rat. Ideal for transcription regulation research.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | SUPT7L Rabbit pAb                                                                                                    |
| Remarks/Alias    | SPT7L; STAF65; SUPT7H; STAF65G; STAF65(gamma); SUPT7L                                                                |
| Species          | Human                                                                                                                |
| Gene ID (Human)  | 9913                                                                                                                 |
| Gene ID          | 9913                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 210-400 of human SUPT7L (NP_001269658) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                    |
| SWISS            | O94864                                                                                                               |
| Protein Weight   | 46kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |