

Anti-Supv3L1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Iyb8

Item Number: CM0024873



Introduction

Rabbit polyclonal antibody targeting SUPV3L1 (SUV3), a mitochondrial RNA helicase. Suitable for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Use in mitochondrial degradosome and RNA metabolism research.

[Learn More](#)

Parameter	Details
Product Name	SUPV3L1 Rabbit pAb
Remarks/Alias	SUV3; SUPV3L1
Species	Human
Gene ID (Human)	6832
Gene ID	6832
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 577-786 of human SUPV3L1 (NP_003162.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IYB8
Protein Weight	88kDa
Shipping	Ice bag