

Anti-Surf1 Rabbit Polyclonal Antibody For Wb And Elisa - Q15526

Item Number: CM0026116

Introduction

Validated for WB and ELISA, this rabbit polyclonal antibody targets human SURF1, a mitochondrial inner membrane protein essential for cytochrome c oxidase assembly. Cross-reacts with mouse and rat, ideal for mitochondrial disease research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	SURF1 Rabbit pAb
Remarks/Alias	SHY1; CMT4K; MC4DN1; SURF1
Species	Human
GeneID (Human)	6834
GeneID	6834
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-275 of human SURF1 (NP_003163.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15526
Protein Weight	33kDa
Shipping	Ice bag