


Anti-Survivin Rabbit Polyclonal Antibody For Wb, Elisa - O15392

Item Number: CM0021626



Introduction

Rabbit polyclonal antibody targeting Survivin (BIRC5/API4), recommended for Western blot and ELISA. Reacts with human, mouse, and rat samples. Ideal for apoptosis and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	Survivin Rabbit pAb
Remarks/Alias	API4; EPR-1; Survivin
Species	Human
Gene ID (Human)	332
Gene ID	332
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human Survivin (NP_001159.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15392
Protein Weight	16 kDa
Shipping	Ice bag