

Anti-Survivin Rabbit Polyclonal Antibody For Wb, If/Icc - O15392

Item Number: CM0023906



Introduction

Rabbit polyclonal antibody against human Survivin (BIRC5, O15392) validated for WB and IF/ICC. Cross-reacts with human and rat. Ideal for apoptosis and cell division research.

[Learn More](#)

Parameter	Value
Product Name	Survivin Rabbit pAb
Remarks/Alias	API4, EPR-1, Survivin(Ab-117)
Species	Human
Gene ID	332 (Human)
Immunogen	Synthetic peptide corresponding to amino acids 1-142 of human Survivin (NP_001159.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC
Cross Reactivity	Human, Rat
SWISS	O15392
Protein Weight	16kDa
Shipping	Ice bag