

Anti-Survivin Rabbit Polymab® Antibody For Wb And Elisa - O15392

Item Number: CM0006506



Introduction

Rabbit monoclonal antibody targeting human Survivin (BIRC5, O15392) validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for apoptosis, cell cycle, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Survivin Rabbit PolymAb®
Remarks/Alias	API4; EPR-1
Species	Human
GeneID (Human)	332
GeneID	332
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human Survivin (NP_001159.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15392
Protein Weight	8-18 kDa
Shipping	Ice bag