

Anti-Suv39H1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43463

Item Number: CM0029849



Introduction

Rabbit polyclonal antibody against human SUV39H1, a histone H3K9 methyltransferase critical for heterochromatin formation and transcriptional repression. Validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	SUV39H1 Rabbit pAb
Remarks/Alias	MG44; KMT1A; SUV39H; H3-K9-HMTase 1; SUV39H1
Species	Human
GeneID (Human)	6839
GeneID	6839
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human SUV39H1 (NP_003164.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43463
Protein Weight	48kDa
Shipping	Ice bag