

Anti-Suv39H2 Polyclonal Antibody For Wb, Elisa, Chip - Q9H5I1

Item Number: CM0024877



Introduction

High-quality rabbit polyclonal antibody targeting human SUV39H2 (KMT1B), validated for WB, ELISA, and ChIP. Ideal for epigenetics research on histone methylation, heterochromatin, and gene silencing.

[Learn More](#)

Parameter	Value
Product Name	SUV39H2 Rabbit pAb
Remarks/Alias	KMT1B; SUV39H2
Species	Human
Gene ID (Human)	79723
Gene ID	79723
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 141-350 of human SUV39H2 (NP_001180354.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H5I1
Protein Weight	47kDa
Shipping	Ice bag