


Anti-Sv2A Rabbit Polyclonal Antibody For Wb And Elisa - Q7L0J3

Item Number: CM0023456



Introduction

Rabbit polyclonal antibody against human SV2A, validated for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Recombinant immunogen targeting amino acids 356-447 (NP_055664.3). Ideal for synaptic vesicle research.

[Learn More](#)

Parameter	Value
Product Name	SV2A Rabbit pAb
Remarks/Alias	SV2, SLC22B1, SV2A
Species	Human
Gene ID	9900 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 356-447 of human SV2A (NP_055664.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q7L0J3
Protein Weight	83kDa
Shipping	Ice bag