

Anti-Svop Rabbit Polyclonal Antibody For Wb, Elisa - Q8N4V2

Item Number: CM0016769

Introduction

High-quality rabbit polyclonal antibody targeting human SVOP (Synaptic vesicle 2-related protein). Validated for Western Blot (WB) and ELISA applications. Ideal for neuroscience research and IVD assay development. Provided by CamelBio.

[Learn More](#)

Parameter	Value
Product Name	SVOP Rabbit pAb
Remarks/Alias	SLC22B4; SVOP
Species	Human
Gene ID	55530
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human SVOP (NP_061181.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8N4V2
Protein Weight	61kDa
Shipping	Ice bag