

Anti-Swap70 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Uh65

Item Number: CM0012984



Introduction

Rabbit polyclonal antibody against human SWAP70, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Target protein weight ~69 kDa. UniProt: Q9UH65.

[Learn More](#)

Parameter	Value
Product Name	SWAP70 Rabbit pAb
Remarks/Alias	HSPC321; SWAP-70; SWAP70
Species	Human
GeneID (Human)	23075
GeneID	23075
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 316-585 of human SWAP70 (NP_055870.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UH65
Protein Weight	69kDa
Shipping	Ice bag