

Anti-Syce1 Polyclonal Antibody For Wb, Elisa - Q8N0S2

Item Number: CM0027499

Introduction

IVD-grade rabbit polyclonal antibody targeting human SYCE1 (CT76) for WB and ELISA. Designed with conserved immunogen, cross-reacts with mouse and rat. Ideal for meiosis and synaptonemal complex studies.

[Learn More](#)

Parameter	Value
Product Name	SYCE1 Rabbit pAb
Remarks/Alias	CT76; POF12; SPGF15; C10orf94; SYCE1
Species	Human
GeneID	93426
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 132-351 of human SYCE1 (NP_001137236.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N0S2
Protein Weight	40kDa
Shipping	Ice bag