

Anti-Sycp3 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Izu3

Item Number: CM0027490



Introduction

Rabbit polyclonal antibody targeting human SYCP3, a testis-specific synaptonemal complex protein essential for meiosis and male fertility. Validated for WB and ELISA, with mouse and rat cross-reactivity. Ideal for reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	SYCP3 Rabbit pAb
Remarks/Alias	COR1; SCP3; SPGF4; RPRGL4; SYCP3
Species	Human
Gene ID (Human)	50511
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-236 of human SYCP3 (NP_710161.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8IZU3
Protein Weight	28kDa
Shipping	Ice bag