

Anti-Syk Monoclonal Antibody For Wb, Ip, If, Fc, Elisa - P43405

Item Number: CM0009235



Introduction

Rabbit monoclonal antibody targeting human Syk, validated for WB, IP, IF, FC, and ELISA. Cross-reacts with mouse and rat. Ideal for BCR signaling, immunoreceptor studies, and hematologic research.

[Learn More](#)

Parameter	Value
Product Name	Syk Rabbit mAb
Remarks/Alias	IMD82; p72-Syk
Species	Human
GeneID (Human)	6850
GeneID	6850
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 251-350 of human Syk (NP_003168.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, FC (intra), ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43405
Protein Weight	72kDa/70kDa
Shipping	Ice bag