

Anti-Syk Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P43405

Item Number: CM0020864



Introduction

Rabbit polyclonal antibody against human SYK (p72-Syk), validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for B-cell receptor signaling and immune research.

[Learn More](#)

Parameter	Value
Product Name	Syk Rabbit pAb
Remarks/Alias	IMD82; p72-Syk; Syk
Species	Human
GeneID	6850 (Human)
Immunogen	Synthetic peptide within amino acids 250-350 of human Syk (NP_001128524.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43405
Protein Weight	72kDa
Shipping	Ice bag