

Anti-Sympk Rabbit Polyclonal Antibody For Wb/Elisa - Q92797

Item Number: CM0029786



Introduction

High-quality rabbit polyclonal antibody targeting human SYMPK (Symplekin). Validated for WB and ELISA, cross-reacts with human and mouse. Ideal for studying histone mRNA processing, tight junction biology, and pre-mRNA polyadenylation.

[Learn More](#)

Parameter	Value
Product Name	SYMPK Rabbit pAb
Remarks/Alias	SPK; SYM; Pta1; SYMPK
Species	Human
GeneID	8189
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-530 of human SYMPK (NP_004810.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92797
Protein Weight	141kDa
Shipping	Ice bag