

Anti-Syn1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P17600

Item Number: CM0028727

Introduction

SYN1 Rabbit Polyclonal Antibody validated for WB, IF-P, ELISA. Reacts with human, mouse, rat SYN1. Ideal for synaptic vesicle trafficking and neurotransmitter release research.

[Learn More](#)

Parameter	Value
Product Name	SYN1 Rabbit pAb
Remarks/Alias	SYN1; EPILX; MRX50; SYN1a; SYN1b; EPILX1; SYN1
Species	Human
GeneID (Human)	6853
GeneID	6853
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human SYN1 (NP_008881.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17600
Protein Weight	74kDa
Shipping	Ice bag