

Anti-Synapsin-1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P17600

Item Number: CM0004477

Introduction

Anti-Synapsin-1 rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA in human, mouse, rat. Targets synaptic vesicle phosphoprotein critical for neurotransmitter release and synaptic plasticity.

[Learn More](#)

Parameter	Value
Product Name	Synapsin-1 Rabbit mAb
Remarks/Alias	SYNI; EPILX; MRX50; SYN1a; SYN1b; EPILX1; Synapsin-1
Species	Human
GeneID (Human)	6853
GeneID	6853
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 401-500 of human Synapsin-1 (NP_008881.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17600
Protein Weight	74kDa
Shipping	Ice bag