

Anti-Synapsin II Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q92777

Item Number: CM0002183

Introduction

Anti-Synapsin II rabbit monoclonal antibody validated for Western blot, immunofluorescence, and ELISA. Reacts with human, mouse, and rat. For synapse and neurotransmitter release research.

[Learn More](#)

Parameter	Value
Product Name	Synapsin II Rabbit mAb
Remarks/Alias	SYNII; Synapsin II
Species	Human
GeneID (Human)	6854
GeneID	6854
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human Synapsin II (Q92777).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92777
Protein Weight	63kDa
Shipping	Ice bag