

Anti-Synaptophysin Polyclonal Antibody For Wb, Ip, If, Ihc - P08247

Item Number: CM0027282

Introduction

Rabbit polyclonal antibody against human Synaptophysin (34kDa) for WB, IP, IF/ICC, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for synaptic vesicle research and neurological disorder studies.

[Learn More](#)

Parameter	Value
Product Name	Synaptophysin Rabbit pAb
Remarks/Alias	MRX96; MRXSYP; XLID96; Synaptophysin
Species	Human
GeneID (Human)	6855
GeneID	6855
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 222-313 of human Synaptophysin (NP_003170.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08247
Protein Weight	34kDa
Shipping	Ice bag