

Anti-Synaptophysin Rabbit Monoclonal Antibody For Wb, If, Elisa - P08247

Item Number: CM0002692

Introduction

High-quality rabbit monoclonal antibody targeting human Synaptophysin. Validated for WB, IF, and ELISA. Cross-reacts with mouse and rat. An essential tool for studying synaptic vesicle biology and plasticity.

[Learn More](#)

Parameter	Value
Product Name	Synaptophysin Rabbit mAb
Remarks/Alias	MRX96; MRXSYP; XLID96; Synaptophysin
Species	Human
GeneID (Human)	6855
GeneID	6855
Immunogen	Synthetic peptide corresponding to amino acids 200-300 of human Synaptophysin (NP_003170.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P08247
Protein Weight	34 kDa
Shipping	Ice bag