

Anti-Synaptophysin Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P08247

Item Number: CM0008990

Introduction

Synaptophysin Rabbit Monoclonal Antibody (Human, Mouse, Rat) validated for WB, IF-P, IHC-P, ELISA. Ideal for studying synaptic vesicle biology, plasticity, and neuropathology. This antibody targets the presynaptic marker synaptophysin, enabling reliable detection in brain tissue and cultured neurons.

[Learn More](#)

Parameter	Value
Product Name	Synaptophysin Rabbit mAb
Remarks/Alias	MRX96; MRXSYP; XLID96
Species	Human
GeneID (Human)	6855
GeneID	6855
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08247
Protein Weight	34kDa
Shipping	Ice bag