

Anti-Syndecan 3 Rabbit Monoclonal Antibody For Wb, Elisa - O75056

Item Number: CM0002351

Introduction

High-quality rabbit monoclonal antibody against human Syndecan-3 (SDC3, O75056) validated for Western blotting and ELISA. Shows cross-reactivity with mouse and rat. Suitable for studying cell surface proteoglycan function and actin cytoskeleton organization.

[Learn More](#)

Parameter	Value
Product Name	Syndecan 3 Rabbit mAb
Remarks/Alias	SDCN; SYND3; Syndecan 3
Species	Human
GeneID	9672
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Syndecan 3 (O75056).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75056
Protein Weight	45kDa
Shipping	Ice bag