

Anti-Syne3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Zmz3

Item Number: CM0018206



Introduction

Rabbit polyclonal antibody against human SYNE3 (Nesprin-3), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Ideal for LINC complex and nuclear positioning studies.

[Learn More](#)

Parameter	Value
Product Name	SYNE3 Rabbit pAb
Remarks/Alias	KASH3; NET53; Nesp3; C14orf49; C14orf139; LINC00341; NCRNA00341; SYNE3
Species	Human
GeneID (Human)	161176
GeneID	161176
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 740-920 of human SYNE3 (NP_689805.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6ZMZ3
Protein Weight	112kDa
Shipping	Ice bag