

Anti-Syngap Rabbit Monoclonal Antibody For Wb And Elisa - Q96Pv0

Item Number: CM0002694

Introduction

Rabbit monoclonal antibody targeting human SynGAP (148 kDa), validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for research on synaptic signaling, plasticity, and autism spectrum disorders.

[Learn More](#)

Parameter	Value
Product Name	SynGAP Rabbit mAb
Remarks/Alias	MRD5; RASA1; RASA5; SYNGAP; SynGAP
Species	Human
Gene ID (Human)	8831
Gene ID	8831
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96PV0
Protein Weight	148kDa
Shipping	Ice bag