


## Anti-Synj1 Polyclonal Antibody For Wb And Elisa - O43426

Item Number: CM0019919



### Introduction

Rabbit polyclonal antibody targeting human SYNJ1 (Synaptojanin-1), validated for Western blot (WB) and ELISA. Detects human and rat SYNJ1 (~173 kDa), a phosphoinositide phosphatase implicated in synaptic endocytosis and Parkinson disease.

[Learn More](#)

| Parameter        | Value                                                                                                 |
|------------------|-------------------------------------------------------------------------------------------------------|
| Product Name     | SYNJ1 Rabbit pAb                                                                                      |
| Remarks/Alias    | DEE53, EIEE53, INPP5G, PARK20, SYNJ1                                                                  |
| Species          | Human                                                                                                 |
| GeneID (Human)   | 8867                                                                                                  |
| GeneID           | 8867                                                                                                  |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 100-200 of human SYNJ1 (NP_003886.3) |
| Source           | Rabbit                                                                                                |
| Category         | Polyclonal Antibodies                                                                                 |
| Application      | WB, ELISA                                                                                             |
| Cross Reactivity | Human, Rat                                                                                            |
| SWISS            | O43426                                                                                                |
| Protein Weight   | 173kDa                                                                                                |
| Shipping         | Ice bag                                                                                               |