

Anti-Synpo Rabbit Polyclonal Antibody For Wb, Elisa - Q8N3V7

Item Number: CM0018327

Introduction

Rabbit polyclonal antibody against human Synaptopodin (SYNPO), validated for WB and ELISA. Detects human, mouse, and rat targets. Useful for studies of actin dynamics in dendritic spines and podocyte foot processes.

[Learn More](#)

Parameter	Value
Product Name	SYNPO Rabbit pAb
Remarks/Alias	SYNPO
Species	Human
Gene ID (Human)	11346
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 820-920 of human SYNPO (NP_001159680.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N3V7
Protein Weight	99kDa
Shipping	Ice bag