

Anti-Synpo2L Rabbit Polyclonal Antibody For Wb, Elisa - Q9H987

Item Number: CM0030031

Introduction

Rabbit anti-SYNPO2L polyclonal antibody validated for WB and ELISA. Detects human and mouse SYNPO2L (~102 kDa), aiding actin-based shape modulation and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	SYNPO2L Rabbit pAb
Remarks/Alias	SYNPO2L
Species	Human
GeneID (Human)	79933
GeneID	79933
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 700-800 of human SYNPO2L (NP_001107605.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9H987
Protein Weight	102kDa
Shipping	Ice bag