

Anti-Syntaphilin (Snph) Rabbit Polyclonal Antibody For Wb, Elisa - O15079

Item Number: CM0013751

Introduction

Rabbit polyclonal antibody targeting human SNPH (Syntaphilin). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for studying SNARE complex inhibition and synaptic function.

[Learn More](#)

Parameter	Value
Product Name	SNPH Rabbit pAb
Remarks/Alias	SNPH; bA314N13.5
Species	Human
Gene ID	9751
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 300-425 of human SNPH (NP_055538.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O15079
Protein Weight	54kDa
Shipping	Ice bag