

Anti-Syntaxin 16 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O14662

Item Number: CM0008967



Introduction

Anti-Syntaxin 16 rabbit monoclonal antibody for WB, IHC-P, ELISA. Validated in human, mouse, rat. Targets a SNARE protein involved in late endosome-to-trans-Golgi vesicular transport. Ideal for membrane trafficking research.

[Learn More](#)

Parameter	Description
Product Name	Syntaxin 16 Rabbit mAb
Remarks/Alias	SYN16; SYN-16; Syntaxin 16
Species	Human
Gene ID	8675
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Syntaxin 16 (O14662)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14662
Protein Weight	37kDa
Shipping	Ice bag