

Anti-Syntaxin 3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13277

Item Number: CM0011852

Introduction

Rabbit monoclonal antibody against human Syntaxin 3 (Q13277), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for studies of membrane trafficking and neurotransmitter transport.

[Learn More](#)

Parameter	Value
Product Name	Syntaxin 3 Rabbit mAb
Remarks/Alias	MVID2, STX3A, DIAR12, RDMVID, Syntaxin 3
Species	Human
GeneID (Human)	6809
GeneID	6809
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human Syntaxin 3 (Q13277)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13277
Protein Weight	33kDa
Shipping	Ice bag