

Anti-Syntaxin (Stx1A/Stx1B/Stx2/Stx3) Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - Q16623 / P61266 / P32856 / Q13277

Item Number: CM0003848

Introduction

Syntaxin rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA. Recognizes human, mouse, and rat syntaxin-1A/1B/2/3. Suitable for SNARE-mediated exocytosis and neuronal function research.

[Learn More](#)

Parameter	Value
Product Name	Syntaxin Rabbit mAb
Remarks/Alias	STX1A/STX1B/STX2/STX3; Syntaxin
Species	Human
GeneID (Human)	6804/112755/2054/6809
GeneID	6804/112755/2054/6809
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Syntaxin (P61266)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16623/P61266/P32856/Q13277
Protein Weight	33 kDa
Shipping	Ice bag