


Anti-Syntenin-1 (Mda-9) Polyclonal Antibody For Wb, If/Icc, Elisa - O00560

Item Number: CM0027275



Introduction

Rabbit polyclonal antibody against human Syntenin-1 (MDA-9) for WB, IF/ICC, and ELISA. Ideal for studying exosome biogenesis, TGF-beta signaling, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	Syntenin Rabbit pAb
Remarks/Alias	ST1, MDA9, SYCL, MDA-9, TACIP18, Syntenin
Species	Human
GeneID (human)	6386
GeneID	6386
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human Syntenin (NP_001007068.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O00560
Protein Weight	32kDa
Shipping	Ice bag