


Anti-Syntenin-1 (Syntenin) Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00560

Item Number: CM0011479



Introduction

Syntenin Rabbit mAb (MDA-9/ST1) for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. High-specificity monoclonal antibody. Explore exosome biogenesis, cancer metastasis, and TGF-beta signaling. | 32kDa | O00560

[Learn More](#)

Parameter	Value
Product Name	Syntenin Rabbit mAb
Remarks/Alias	ST1; MDA9; SYCL; MDA-9; TACIP18; Syntenin
Species	Human
GeneID	6386
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Syntenin (O00560).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00560
Protein Weight	32kDa
Shipping	Ice bag