

Anti-Syntrophin Alpha 1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13424

Item Number: CM0023626

Introduction

Rabbit polyclonal antibody against human Syntrophin alpha 1 (SNTA1). Validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying dystrophin complex and neuromuscular synapse.

[Learn More](#)

Parameter	Value
Product Name	Syntrophin alpha 1 Rabbit pAb
Remarks/Alias	SNT1; LQT12; TACIP1; dj1187j4.5; Syntrophin alpha 1
Species	Human
Gene ID (Human)	6640
Gene ID	6640
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human Syntrophin alpha 1 (NP_003089.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13424
Protein Weight	54kDa
Shipping	Ice bag