

Anti-Sypl1 Polyclonal Antibody For If/Icc, Elisa - Q16563

Item Number: CM0021271

Introduction

Anti-SYPL1 rabbit polyclonal antibody validated for IF/ICC and ELISA applications. Reacts with human synaptophysin-like protein 1 (pantophysin), a membrane protein involved in cytoplasmic vesicle traffic. Ideal for cell biology research.

[Learn More](#)

Parameter	Value
Product Name	SYPL1 Rabbit pAb
Remarks/Alias	SYPL; H-SP1; SYPL1
Species	Human
Gene ID (Human)	6856
Gene ID	6856
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-120 of human SYPL1 (NP_006745.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q16563
Protein Weight	29kDa
Shipping	Ice bag