

Anti-Syt1 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P21579

Item Number: CM0008059



Introduction

High-performance rabbit monoclonal antibody targeting human SYT1 (synaptotagmin-1). Validated for WB, IF/ICC, IHC-P, ELISA, with cross-reactivity to mouse and rat. Ideal for synaptic vesicle research.

[Learn More](#)

Parameter	Value
Product Name	SYT1 Rabbit mAb
Remarks/Alias	P65; SYT; BAGOS; SVP65
Species	Human
Gene ID (Human)	6857
Gene ID	6857
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 141-422 of human SYT1 (NP_005630.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21579
Protein Weight	48kDa
Shipping	Ice bag