

Anti-Syt11 Rabbit Polyclonal Antibody For Wb, If, Elisa - Q9Bt88

Item Number: CM0029412

Introduction

Rabbit polyclonal antibody against human Synaptotagmin-11 (SYT11), validated for WB, IF/ICC, IF-P, ELISA. Detects endogenous SYT11 in human, mouse, and rat samples, with a predicted molecular weight of 48 kDa.

[Learn More](#)

Parameter	Value
Product Name	SYT11 Rabbit pAb
Remarks/Alias	SYT12; sytXI; SYT11
Species	Human
GeneID (Human)	23208
GeneID	23208
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-200 of human SYT11 (NP_689493.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BT88
Protein Weight	48kDa
Shipping	Ice bag