


Anti-Syt4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9H2B2

Item Number: CM0027808



Introduction

Rabbit polyclonal antibody targeting human Synaptotagmin-4 (SYT4), validated for WB, IHC-P, IF/ICC, and ELISA. Shows cross-reactivity with mouse and rat. Suitable for neuronal dense core vesicle trafficking and melanocyte dendrite formation research.

[Learn More](#)

Parameter	Value
Product Name	SYT4 Rabbit pAb
Remarks/Alias	HsT1192; SYT4
Species	Human
GeneID (Human)	6860
GeneID	6860
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-190 of human SYT4 (NP_065834.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H2B2
Protein Weight	48kDa
Shipping	Ice bag