

Anti-Syt6 Rabbit Polyclonal Antibody For Wb, Elisa - Q5T7P8

Item Number: CM0020484

Introduction

Anti-SYT6 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat synaptotagmin-6, a calcium sensor in vesicular trafficking and exocytosis. Ideal for neuroscience and secretion research.

[Learn More](#)

Parameter	Value
Product Name	SYT6 Rabbit pAb
Remarks/Alias	sytVI; SYT6
Species	Human
GeneID (Human)	148281
GeneID	148281
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human SYT6 (NP_001257734.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5T7P8
Protein Weight	57kDa
Shipping	Ice bag