

Anti-Syt7 Rabbit Polyclonal Antibody For Wb And Elisa - O43581

Item Number: CM0014459



Introduction

Rabbit polyclonal antibody against human Synaptotagmin-7 (SYT7), validated for Western blot (WB) and ELISA. SYT7 is a calcium sensor essential for lysosomal exocytosis, membrane repair, and synaptic facilitation.

[Learn More](#)

Parameter	Value
Product Name	SYT7 Rabbit pAb
Remarks/Alias	IPCA7; IPCA-7; SYTVII; PCANAP7; SYT-VII; SYT7
Species	Human
Gene ID (Human)	9066
Gene ID	9066
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-170 of human SYT7 (NP_004191.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43581
Protein Weight	46kDa
Shipping	Ice bag