


Anti-T1R3/Tas1R3 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q7Rtx0

Item Number: CM0007342



Introduction

High-performance rabbit monoclonal antibody targeting T1R3/TAS1R3, a key taste receptor subunit for sweet and umami perception. Validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	T1R3/TAS1R3 Rabbit mAb
Remarks/Alias	T1R3
Species	Human
GeneID	83756
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-570 of human T1R3/TAS1R3 (NP_689414.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7RTX0
Protein Weight	93kDa
Shipping	Ice bag