


Anti-Tab2 (Map3K7Ip2) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nyj8

Item Number: CM0027433



Introduction

Explore anti-TAB2 rabbit polyclonal antibody for WB and ELISA applications. Targets human, mouse, and rat TAB2 (UniProt Q9NYJ8), a key adapter in NF- κ B and JNK signaling.

[Learn More](#)

Parameter	Value
Product Name	TAB2 Rabbit pAb
Remarks/Alias	CHTD2; TAB-2; MAP3K7IP2; TAB2
Species	Human
Gene ID (Human)	23118
Gene ID	23118
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-620 of human TAB2 (NP_055908.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYJ8
Protein Weight	76kDa
Shipping	Ice bag