

Anti-Tab3 Monoclonal Antibody For Wb, Elisa - Q8N5C8

Item Number: CM0003650

Introduction

High-specificity rabbit monoclonal antibody targeting TAB3 (NAP1/MAP3K7IP3), validated for Western blot and ELISA in human and mouse samples. Ideal for studying NF-kB and JNK signaling in cancer and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	TAB3 Rabbit mAb
Remarks/Alias	NAP1; MAP3K7IP3; TAB3
Species	Human
GeneID (Human)	257397
GeneID	257397
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human TAB3 (Q8N5C8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N5C8
Protein Weight	79kDa
Shipping	Ice bag