

Anti-Tac1 Rabbit Polyclonal Antibody For Wb, Elisa - P20366

Item Number: CM0019066



Introduction

Anti-TAC1 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human and rat TAC1. Recombinant immunogen corresponding to amino acids 1-129. Ideal for neuroscience and inflammation research involving substance P and tachykinin signaling.

[Learn More](#)

Parameter	Value
Product Name	TAC1 Rabbit pAb
Remarks/Alias	NK2; NPK; NKNA; TAC2; Hs.2563; TAC1
Species	Human
Gene ID (Human)	6863
Gene ID	6863
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-129 of human TAC1 (NP_003173.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human,Rat
SWISS	P20366
Protein Weight	15 kDa
Shipping	Ice bag